



USB
Charger
Port



ELEGANCE COLLECTION

Now with
FIRST
TECHNOLOGY

UP 50%
Faster
LIFT
SPEED



Scan QR code
for video!

Available on the LC-400 Series

The Elegance Collection from Pride® Lift Chairs is designed to be the ultimate in lift chair comfort, style and performance. The state-of-the-art lift system used in Pride Elegance lift chairs allows for a full recline, making it the ideal solution for napping and relaxation.

AMERICA'S #1
LIFT CHAIR BRAND

Pride®
Lift Chairs

1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)
pridemobility.com • prideliftchairs.com

Pride®
Mobility Products Corp.
— LIVE YOUR BEST —

400 SERIES




SPECIFICATIONS

	Medium	Medium	Medium	Small	Medium	Medium Wide	Large Tall	Medium
	LC-421 ^{2,4}	LC-450 ^{2,4}	LC-465	LC-470S ⁴	LC-470M ^{2,4}	LC-470W ^{2,4}	LC-470LT ^{2,4}	LC-485 ^{2,4}
	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 6'1" - 6'6"	3-Position Chaise Lounger 5'5" - 6'0"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	20"	20"	20"	19"	20"	20"	21"	20"
Seat Depth	20"	20"	20"	19"	20"	20"	22"	20"
Seat Width	20"	20"	20"	19"	20"	24"	22"	20"
Top of Back to Seat	24.5"	27"	26"	24.5"	26.5"	26.5"	31"	26"
Back Style	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Pillow Back	Blown Fiber Pillow Back	Blown Fiber Waterfall Back	Blown Fiber Waterfall Back	Blown Fiber Waterfall Back	Blown Fiber Biscuit Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Total Comfort Seating	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Distance From Wall	19"	17"	17"	16.5"	19"	19"	24"	18"
Overall Chair Width	38"	38"	36"	35"	36"	40"	41"	36"
Chair Weight	119 lbs.	119 lbs.	113 lbs.	110 lbs.	113 lbs.	117 lbs.	114 lbs.	120 lbs.
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated

DURASOFT FABRICS



CLOUD 9 FABRICS



CLOUD 9 OPTIONAL FABRICS



PADDED-SUEDE



OPTIONAL FABRICS



LEXIS URETHANE



ULTRAFABRICS® ULTRALEATHER™ COLORS



CRYPTON® SUPER FABRICS

SMART SUEDE³



SMART FABRICS³



OPTIONAL MESA FABRICS



BENEFITS

- Industry-first hand control with USB charger port (not available with Heat and Massage)
- Standard head & arm covers
- Standard large dual pockets
- Removable back
- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (except LC-421)
- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Fire retardant (California specifications; bulletin 117-2013)
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 9' cord
- Lifetime warranty on steel lift frame
- 7-year prorated warranty

OPTIONS*

- Deluxe Heat & Massage
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levers
- Footrest extension available (except LC-465 & LC-470S)
- First Up technology

Additional lead time may
apply for some options.

*Additional charges apply for all options.

1. Batteries not included.
2. Available with footrest extension.
3. Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.
4. Available with Deluxe Heat & Massage.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



Look for the Quick Ship
Indicator for the shortest
lead times available.