



TECHNICAL SPECIFICATIONS

Standard Seat Pad Size: 17.5" / 445 mm wide x 15" / 381 mm deep

Larger Seat Pad Option: 21.5" / 546 mm wide x 16" / 406 mm deep

Standard Footrest Size: 10" / 254 mm long x 9" / 229 mm wide

Larger Footrest Option: 12" / 305 mm long x 11" / 279 mm wide



High-quality American craftsmanship is the hallmark of custom-made Bruno Elite curved rail stairlifts. Designed to the exact specifications of your staircase, the Elite melds into every curve. The plush chair properly reflects the beauty of your home. Let Bruno's Elite safely connect you to all levels of your home again.

- Weight capacity: 400 lb/ 181 kg
- Plush, generous-size seat with multiple height adjustments
- Adjustable width arm rests to 21.75" / 552 mm between arms
- Offset swivel seat makes getting on/off safe and easy
- Under-armrest rocker switch for easy operation
- Painted rail with covered gear rack gives a sleek appearance
- Flip up arms, seat, and footrest create plenty of extra space on steps
- Two wireless call/sends allow unit to be controlled remotely
- Two 12-v batteries ensure dependable performance even in a power outage
- Soft start/stop for a smooth ride from start to finish
- Adjustable footrest height for maximum comfort
- Onboard audio/visual diagnostics for easy service
- Obstruction safety sensors, seat belt ensure safety
- Gold warranty: five year major components; two-year parts

See reverse side for optional features. Learn more at www.bruno.com

ELITE CURVED




ELITE CURVED OPTIONS

Standard tan vinyl with five optional upholstery selections. Or, supply your own fabric and Bruno will custom upholster the chair for you.

• Wider seat

• Wider footrest



Top or bottom park position extends rail away from the staircase.



Power folding footrest automatically flips up/down when seat is raised/lowered.



Power swivel seat for effortless exit. Controlled on chair arm or wireless call/send.



Mid-park and charge station for staircases with middle landings.



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2015

MKT-SL-3