

# Guardian Wheelchairs

*As the world's leading manufacturer of custom mobility products, Sunrise Medical knows wheelchairs. We also know that for basic patient transportation needs your requirements are different. You need wheelchairs that have been designed with durability in mind to provide years of trouble-free service at a cost that maximizes the return on your investment.*

*You also require choices. Choices in style, size, features and configurations to accommodate all of your customers' needs, as well as a choice of suppliers.*

*Six Easy Care™ models, covering the spectrum from basic to high-strength lightweight, along with the Guardian Escort®, the new standard for the 21st century, establish the Guardian family of wheelchairs as the best combination of quality, value and selection.*



## Guardian Escort®


*The objective was to create a standard wheelchair combining quality, function and value to elevate the industry standard for the 21st century. We're proud to say the Guardian Escort exceeds these expectations. Designed with technologically advanced materials, the results are a product with an elegantly simple design and rugged durability.*

### Standard Features

- Innovative design using technologically advanced materials to provide a scratch-resistant, rust-resistant and extremely durable chair
- Distinctive style and attractive silver texture finish is easy to maintain
- Armrests feature a 3-point anchor and toggle release for secure attachment and easy removal
- Spring-assisted wheel lock design provides improved user safety during transfers
- Adjustable axle bracket improves maneuverability and stability

### Options for Versatility

- Bi-directional footrests swing in and out for easier transfers and better access in small spaces
- Unique elevating legrest design eliminates traditional push rod to deliver long-term performance
- Full length, traditional or standard desk length arms with armpads made of advanced composite materials provide the comfort of padding and the ease of plastic
- Transport chair configuration with 12" rear wheels and fold-down space-saving backrest
- Padded nylon or easy-to-clean vinyl upholstery provides comfort and convenience



**Swing In/Out Footrests**



**Adjustable Axle**

# Specifications

All specifications are subject to change without notice.

	Easy Care 1000 Transport	Easy Care 1000	Easy Care 2000	Easy Care 2000HD	Easy Care 3000	Easy Care 4000	Easy Care 4000 Recliner	Escort Standard
Seat Width	17"/19"	16"/18"	16"/18"	20"/22"/24"	16"/18"	16"/18"/20"	18"	16"/18"/20"
Seat Depth	16"	16"	16"	16"/18"/18"	16"	16"/18"	17"	16"
Overall Width								
Fixed Arm	21.5"/23.5"	22.5"/24.5"	22.5"/24.5"	N/A	N/A	N/A	N/A	N/A
Removable Arm	N/A	N/A	24.5"/26.5"	29.5"/30.5"/32.5"	23"/25"	23"/25"/27"	25"	24.5"/26.5"/28.5"
Overall Depth								
Fixed Arm, Fixed F/R	N/A	41.5"	41.5"	N/A	N/A	N/A	N/A	N/A
Fixed Arm w/o riggings	N/A	N/A	32"	N/A	N/A	N/A	N/A	N/A
Removable w/o riggings	24"	N/A	29"	29"/32"/32"	30.5"	30.5"	34.5"	42"
Seat to Floor Height	19"	19.75"	19.75"	19.75"	19.75"	19.75"	19.75"	20"
Back Style	Fold Down	Fixed	Fixed	Fixed	Fixed	Adjustable	Recline	Fixed
Back Height	17"	16"	16"	16"	16"	16" - 18"	22"	17"
Arm Styles	Fixed Full	Fixed Full	Fixed Full Removable Desk Removable Full	Removable Full Removable Desk	Removable Desk	Flip-Back Desk Flip-Back Full	Removable Desk	Removable Full Traditional Desk Standard Desk
Front Riggings	Swing-away	Fixed	Fixed Swing-away Elevating	Swing-away Elevating	Swing-away Elevating	Swing-away Elevating	Swing-away Elevating	Swing-away Elevating
Rear Axle	Single Position	Single Position	Dual Position	Dual Position (20") Single Posit. (22", 24")	Dual Position	Dual Position	Dual Position	Multiple Axle
Rear Wheels	8x1 in. Rubber	Composite Mag	Composite Mag	Composite Mag	Composite Mag	Composite Mag	Composite Mag	Composite Mag
Hand Rims	N/A	Composite	Composite	Chromed Steel	Composite	Aluminum	Aluminum	Painted Steel
Upholstery	Nylon	Vinyl	Vinyl	Vinyl	Nylon	Nylon	Nylon	Nylon
Product Weight (w/o riggings)	19 lbs. (Aluminum) 26 lbs. (Steel)	41 lbs.	41 lbs.	41 lbs. (20") 57 lbs. (22"/24")	35.5 lbs.	33.5 lbs.		37 lbs.
Weight Capacity	250 lbs.	250 lbs.	250 lbs.	300/450/450 lbs.	250 lbs.	250 lbs.	250 lbs.	250 lbs.
Warranty								
Frame	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	5 Years
Tires	3 Years	3 Years	3 Years	3 Years	3 Years	3 Years	3 Years	N/A
Upholstery	90 Days	90 Days	90 Days	90 Days	90 Days	90 Days	90 Days	N/A
Other Components	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	3 Years
Reimbursement Code <sup>1</sup>	E1031	K0001	K0001	K0006/K0007 K0006	K0003 K0001/K0002	K0004 K0001/K0002/K0003	K0004 & K0028	K0001
Eligibility For No Charge Upgrade								

## Accessories

	Easy Care 1000 Transport	Easy Care 1000	Easy Care 2000	Easy Care 2000HD	Easy Care 3000	Easy Care 4000	Easy Care 4000 Recliner	Escort Standard
Overhead Anti-Theft Device	N/A	GEY-OATH	GEY-OATH	GEY-OATH	GEY-OATH	GEY-OATH	N/A	N/A
Rear Anti Tippers	N/A	GEY-AT	GEY-AT	GEY-AT	GEY-AT	GEY-4AT	Included	GE-AT
Anti-Fold Device	N/A	GEY-ATH	GEY-ATH	GEY-ATH	GEY-ATH	GEY-ATH	GEY-ATH	GE-AF
Oxygen Tank Holder <sup>2</sup>	N/A	GEY-O2	GEY-O2	GEY-O2	GEY-O2	GEY-O2	GEY-O2	GE-O2/GE-O2IV
Telescoping IV Pole	N/A	GEY-IVP	GEY-IVP	GEY-IVP	GEY-IVP	GEY-IVP	GEY-IVP	GE-IV
Wheel Lock Extension Handles	N/A	GEY-EXT	GEY-EXT	GEY-EXT	GEY-EXT	GEY-EXT	GEY-EXT	GE-6E
Amputee Axle Plate	N/A	GEY-AMP	GEY-AMP	GEY-HDAMP	GEY-AMP	GEY-4AMP	N/A	GE-AMPT
Amputee Axle Plate Wheel Lock	N/A	GEY-WDG	GEY-WDG	GEY-WDG	GEY-WDG	GEY-WDG	N/A	N/A
Chart Pocket	GE-CP	GE-CP	Included	Included	GE-CP	GE-CP	GE-CP	GE-CP
Composite Footrests	GEY-FRC	GEY-FRC	GEY-FRC	GEY-FRC	GEY-FRC	GEY-FRC	GEY-FRC	N/A
Composite Footrests w/Heel Loop	N/A	N/A	GEY-FRCHL	GEY-FRCHL	GEY-FRCHL	GEY-FRCHL	GEY-FRCHL	ESQSFR16/18
Composite Elevating Legrest	N/A	N/A	GEY-ELRC	GEY-ELRC	GEY-ELRC	GEY-ELRC	GEY-ELRC	ESQSEL16/18
Low Profile ELR	N/A	N/A	GEY-ELRCLP	N/A	N/A	GEY-ELRCLP	GEY-ELRCLP	N/A
Articulating Legrests	N/A	N/A	GEY-ALR	GEY-ALR	GEY-ALR	GEY-ALR	GEY-ALR	N/A
Ht. Adj Padded Desk Armrest	N/A	N/A	GEY-HADP	GEY-HADP	GEY-HADP	GEY-4HADP	N/A	ESQSSH <sup>3</sup>
Ht. Adj Padded Full Armrest	N/A	N/A	GEY-HAFP	GEY-HAFP	GEY-HAFP	GEY-4HAFP	N/A	N/A
2" Seat Extension	N/A	N/A	GEY-2SREXT* * = 16, 18 or 20	N/A	N/A	GEY-4SREXT* * = 16 or 18	N/A	N/A
Seat Belt (Automotive Style)	GEY-SBA	GEY-SBA	GEY-SBA	GEY-SBA	GEY-SBA	GEY-SBA	GEY-SBA	N/A
Seat Belt (Velcro Style)	GEY-SBV	GEY-SBV	GEY-SBV	GEY-SBV	GEY-SBV	GEY-SBV	GEY-SBV	N/A

1. Meets requirements for HCPCs code

2. Must use with rear anti-tippers

3. Available as upgrade option on chair



# Guardian Standard Wheelchairs



# GUARDIAN<sup>®</sup>

The best combination of value,  
durability and selection.

## Guardian Easy Care™ 1000 Transport Wheelchair

*Compact and lightweight, this economical model offers the most convenience for transporting patients with the assistance of an attendant.*

### Features

- Back folds down making transport and storage easier
- Choice of chrome-plated carbon steel or anodized aluminum frame
- Lightweight, the aluminum models weigh only 19 lbs. (8.6 kg.)—25% lighter than the steel configuration which is 26 lbs. (11.7 kg) (without footrests)
- Padded nylon upholstery is comfortable and durable
- 8" wheels in front and back provide a smooth ride
- Full-length arms and a nylon seat belt are standard

### Models

- |            |                                |
|------------|--------------------------------|
| IT17FAFPS  | 17" seat width, steel frame    |
| ITA17FAFPS | 17" seat width, aluminum frame |
| IT19FAFPS  | 19" seat width, steel frame    |
| ITA19FAFPS | 19" seat width, aluminum frame |



## Guardian Easy Care™ 1000 Wheelchair

*An economical choice for rugged durability, this no-frills model will complement a rental fleet.*

### Features

- Chrome-plated, carbon steel frame is built to last and easy to maintain
- Composite mag wheels are molded with contours to enhance strength and durability
- Sealed bearings are maintenance free
- Fixed full-length armrests and fixed footrests are standard
- Easy-to-clean vinyl upholstery is flame retardant

### Models

- |           |                |
|-----------|----------------|
| 116FAFUFS | 16" seat width |
| 118FAFUFS | 18" seat width |



Dual axle allows seat to be easily lowered to hemi height. Available on 2000, 3000 and 4000 models.



Chart pocket is included on 2000 and 2000HD models.



Quick release flip-back arms make transfers easy. Available on 4000 models only.



Height adjustable back. Standard on 4000 models only.



# Guardian Easy Care™ 2000 Wheelchair

*A good choice for a quality wheelchair with options that offer the flexibility to accommodate a variety of users.*

## Standard Features

- Dual axle allows the seat to be easily lowered to hemi height
- Chrome-plated, carbon steel frame is built to last
- Composite mag wheels are molded with contours to enhance strength and durability
- Full size chart pocket on backrest
- Front rigging mounts are compatible with the industry's standard

## Options for Versatility

- Removable desk-length or full-length or fixed full-length armrests
- Swing-away footrests, elevating legrests or fixed footrests

## Models

216FAFPFS	16" width, fixed padded full arms, fixed footrests
216FAFPS	16" width, fixed padded full arms, swing-away footrests
216RADPS	16" width, removable padded desk arms, swing-away footrests
216RADPE	16" width, removable padded desk arms, elevating legrests
216RADUS	16" width, removable desk arms, swing-away footrests
216RAFPFS	16" width, removable padded full arms, swing-away footrests
216RAFPE	16" width, removable padded full arms, elevating legrests
218FAFPFS	18" width, fixed padded full arms, fixed footrests
218FAFPS	18" width, fixed padded full arms, swing-away footrests
218RADPS	18" width, removable padded desk arms, swing-away footrests
218RADPE	18" width, removable padded desk arms, elevating legrests
218RADUS	18" width, removable desk arms, swing-away footrests
218RAFPFS	18" width, removable padded full arms, swing-away footrests
218RAFPE	18" width, removable padded full arms, elevating legrests



# Guardian Easy Care™ 2000HD Heavy-Duty and Extra Wide Wheelchairs

*Extra wide seat widths and reinforced construction accommodate special needs.*

## Features

- Heavy-duty reinforced chromed steel frame supports weights up to 450 lbs. (204 kg) (22" and 24" seat width only)
- Available in 20" heavy-duty and 22" or 24" extra wide seat widths
- Durable, double-embossed vinyl upholstery
- Dual axle allows the seat to be easily lowered to hemi height (20" seat width only)
- Swing-away footrests or elevating legrests with padded calf
- Oversize chart pocket

## Models

220RADUS	20" width, desk arms, swing-away footrests
220RADPS	20" width, padded desk arms, swing-away footrests
220RADPE	20" width, padded desk arms, elevating legrests
220RAFPFS	20" width, padded full arms, swing-away footrests
220RAFPE	20" width, padded full arms, elevating legrests
2HD22RADPS	22" width, padded desk arms, swing-away footrests
2HD22RADPE	22" width, padded desk arms, elevating legrests
2HD22RAFPFS	22" width, padded full arms, swing-away footrests
2HD22RAFPE	22" width, padded full arms, elevating legrests
2HD24RADPS	24" width, padded desk arms, swing-away footrests
2HD24RADPE	24" width, padded desk arms, elevating legrests
2HD24RAFPFS	24" width, padded full arms, swing-away footrests
2HD24RAFPE	24" width, padded full arms, elevating legrests



# Guardian Easy Care™ 3000 Lightweight Wheelchair

*A value-priced lightweight design to meet today's funding requirements.*

## Features

- Weighs 35 pounds (16 kg) without riggings meeting requirements for K0003 funding
- Scratch-resistant composite handrim
- Chrome-plated, carbon steel frame is durable
- Dual axle allows seat to be easily adjusted to hemi height
- Removable desk arms with padded armrests
- Padded nylon upholstery is comfortable
- Swing-away footrests or elevating legrests are interchangeable with the industry's standard

## Models

316RADPS	16" width, swing-away footrests
316RADPE	16" width, elevating legrests
318RADPS	18" width, swing-away footrests
318RADPE	18" width, elevating legrests
320RADPS	20" width, swing-away footrests
320RADPE	20" width, elevating legrests



# Guardian Easy Care™ 4000 High-Strength Lightweight Wheelchair

*Durability in a lightweight design and unique options make this the most versatile Easy Care model.*

## Features

- Frame design eliminates seat guides and allows for custom backs and accessories
- Flip-back arms with an easy release mechanism
- 33.5 pound (15 kg) weight without riggings meets the needs of those requiring a high-strength lightweight standard wheelchair
- Height adjustable seat back (16"/17"/18") offers greater comfort

## Models

416RADPE	16" width, padded desk arms, elevating legrests
416RADPS	16" width, padded desk arms, swing-away footrests
416RAFPE	16" width, padded full arms, elevating legrests
416RAFPS	16" width, padded full arms, swing-away footrests
418RADPE	18" width, padded desk arms, elevating legrests
418RADPS	18" width, padded desk arms, swing-away footrests
418RAFPE	18" width, padded full arms, elevating legrests
418RAFPS	18" width, padded full arms, swing-away footrests
420RADPE	20" width, padded desk arms, elevating legrests
420RADPS	20" width, padded desk arms, swing-away footrests
420RAFPE	20" width, padded full arms, elevating legrests
420RAFPS	20" width, padded full arms, swing-away footrests
4R16RADPE	16" width, reclining back, padded desk arms, elevating legrests
4R16RADPS	16" width, reclining back, padded desk arms, swing-away footrest
4R18RADPE	18" width, reclining back, padded desk arms, elevating legrests
4R18RADPS	18" width, reclining back, padded desk arms, swing-away footrest
4R20RADPE	20" width, reclining back, padded desk arms, elevating legrests
4R20RADPS	20" width, reclining back, padded desk arms, swing-away footrests

