



Harmar's Value Stair Lift Gets the Job Done

Designed for cost-conscious buyers; the Vantage Stair Lift is a sturdy, comfortable, low maintenance lift that is easy to install. Utilizing a well protected rack & pinion drive system with a speed regulator, users will experience a smooth, consistent drive speed as they ascend or descend. And with an upholstered seat and powder-coated chassis, the stylish stair lift is sure to complement any home decor.

Features & Benefits

- Speed regulator offers a smooth comfortable ride at a constant speed in both directions
- Modern styling with upholstered seat
- Traditional rack & pinion drive system for cost effectiveness
- Quick & easy installation
- Covered gear rack keeps your stairs clean, looks stylish
- Compact, folds to 13.6" to allow others to pass by
- DC drive system can operate several cycles during power outages
- Powerful for stair inclines up to 45°
- Includes 2 wireless remotes
- ETL Listed - 3148125

Specifications	SL400
Capacity	300 lbs.
Speed	19 feet/min.
Swivel	90° at top
Rail Install	5" from wall; mounts on steps
Drive System	Rack & Pinion
Braking	Self-locking Gear Box
Charge Strip	3 ft. at top & bottom
Folded Width	13.6"
Clear Distance Between Arm Rests	19"
Folding Rail	Not Available
Rail Length	16' - 0" Standard 32' - 0" Maximum
Motor Output	14 amps
Warranty	3-Year Parts 1-Year on Batteries
Safety Design	Complies with ASME A18.1, CAN/CSA B44.1, ASME A17.5
Safety Devices	Footrest Obstruction Sensor Swivel Seat Cut-off Switch Final Limits / Seat Belt