

The popular Buzzaround Lite scooter comes standard with a Delta tiller with wraparound handle and a larger, more comfortable seat. The Buzzaround Lite still features our easy "Pop! Pull! Push!" wireless

disassembly for transport. Change the look of your scooter with our three sets of snap-on shroud panels in Red, Blue and Sahara, all while

enjoying its tight indoor maneuverability and smooth outdoor use.

technologies

2

Rise. Roll.

BUZZ around

<u>lite</u><

3 wheel

GB-106

Sahara

YOUR LIFE IN MOTION



LIGHTWEIGHT PORTABLE SCOOTER

3 wheel Specifications

Model	GB106
Medicare Code	K0800
Weight Capacity	250 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range ¹	8 miles
Ground Clearance to center deck	1.75″
Ground Clearance to transaxle/motor	1.5″
Turning Radius	31″
Type Batteries	2-12AH
Colors	Red, Blue, Sahara
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	45 Amp
Dimensions: ³	
Length	36.5″
Width	21.25″
Height (ground to top of back of seat)	30.375" - 32.375"
Ground to Top of Seat	19.25" - 21.25"
Ground to Top of Deck	4.25″
Deck to Top of Seat	15" - 17"
Number of Seat Height Positions	3
Size of Increments	1″
Front Axle to Rear Axle	28.625″
Weight of Unit: (assembled)	96 lbs
Front Half	29 lbs
Rear Half	26 lbs
Seat w/Arms	21 lbs
Batteries	19 lbs/pack
Tires:	Solid
Front	8″
Rear	8″
Rear Anti-Tippers	2″
Standard Seat:	Stadium Style
Back Height - no headrest	16"
Width x Depth	17″ X 16″
Color	Black Vinyl
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years

^{1.} Battery range will vary due to rider weight, drive surface and terrain.

^{2.} Electronics warranty excludes batteries.

^{3.} There is a +1 "/-1" tolerance on all dimensions.



Exclusive interlocking teeth provide the tightest, most dependable battery connection for lasting power.

Wireless, quick disconnect in three easy steps:





Step 2: Pull off the battery pack!



Step 1: Pop off the seat!

Step 3: Push apart the frame!



Quickly disassembles into five, lightweight pieces for easy transport!

Features

- Exclusive interlocking teeth for fully connectorless assembly/ disassembly
- **250 lb. weight capacity**
- Interchangeable color panels included with each unit–FREE!
- Auto grade cut pile carpet
- Tiller mounted basket



New control panel design and Delta tiller standard



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com