

Seat Color Black Vinyl







3 wheel GB117S 4 wheel GB147S

The next generation of the popular Buzzaround scooter series is now available with our comfort spring system! The Buzzaround XLS offers more standard features to make your riding experience better than ever. Features a wraparound Delta tiller, a new LED battery gauge and a new LED angleadjustable headlight. The new Buzzaround XLS is also longer to offer more foot room, with a great turning radius. The Buzzaround XLS disassembles in seconds with just four easy steps making transporting it from one destination to another a breeze!

 \equiv

S

 \pm

 \simeq Д

 \pm

5

 \Box





LIGHTWEIGHT PORTABLE SCOOTER

Item	GB117S	GB147S
Description	XLS	XLS
Medicare Code	N/A	N/A
Weight Capacity	300 lbs.	300 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to ¹	8 miles	8 miles
Under Deck Clearance	3"	3"
Turning Radius	37"	47"
Type Batteries	2-12AH	2-12AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brake	es Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp Dynamic R-Series	50 Amp Dynamic R-Series
Dimensions		
Length	42.5"	41.5"
Width	21.5"	21.5"
Height (ground to top of back of seat)	33"-35"	33"-35"
Ground to Top of Seat	21"-23"	21"-23'
Ground to Top of Deck	4.75"	4.75'
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height P	ositions 3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	33.5"	32"
Weight of Unit: (assembled)³	104 lbs.	112 lbs
Front Half	35 lbs.	43 lbs
Rear Half	28 lbs.	28 lbs
Seat w/Arms	22 lbs.	22 lbs
Batteries	19 lbs. pack	19 lbs. pack
Tires:	Solid	Solic
Front	8"	8'
Rear	8"	8"
Rear Anti-Tip Wheels	2"	2'
Standard Seat	Stadium Style	Stadium Style
Seat Back Height - no headrest	16"	16'
Width x Depth	18" x 16"	18" x 16'
Color	Black Vinyl	Black Viny
Optional Seat:	Stadium Style	Stadium Style
Back Height - no headrest	16"	16′
Width x Depth	20"x 18"	20"x 18'
Color	Black Vinyl	Black Viny
	Didek viriyi	Didek VIIIy
Warranty:		
Warranty: Frame	Lifetime Warranty	Lifetime Warrant
Frame	Lifetime Warranty	Lifetime Warranty
	Lifetime Warranty 2 Year Warranty 2 Year Warranty	Lifetime Warranty 2 Year Warranty 2 Year Warranty

Battery range tested at 200 lbs., and will vary due to rider weight, drive surface and terrain.






Comfort Spring System

Wireless disassembly in 4 easy steps!







Step 2: Pull out the battery pack!

Step 3: Fold down tiller and lock in place!



Step 4: Pull forward on the drivetrain release lever to separate the front and rear frame sections!




Disassembles into 5 pieces!

4 wheel model shown

Features

- 300 lb. weight capacity
- **■** Tiller mounted basket
- Angle-adjustable LED headlight
- Comfort Spring System





^{2.} Electronics warranty excludes batteries.