



Quickie 2

By Sunrise



Description

The Quickie 2 from Sunrise Medical has more features and options than any other folding chair in today's market. There are over 64,000 configurations to customize your chair to your individual needs. It's versatility and well-established track record has made the Quickie 2 the premier choice of manual folding ultralight wheelchairs for over 25 years.

Picture shows chair with cushioned back and laterals.

Features

- Multi-adjustable caster forks allow for the most possible seat-to-floor height choices (from 13" - 22")
- New swing in and swing out detachable footrests for easier transfer accomidation
- Q-Fit Technology for the most quiet and stable armrest and side guard recievers
- Quickie's Tru-Fit Program provides one-time free parts offer to modify the size of the chair within five years of ownership
- Armrest options allow optimal user support, positioning and transfer needs
- 265 lbs. weight capacity with a 350 lbs. capcity heavy duty option
- Free Shipping

Quickie 2 Foldable Ultra Lightweight Wheelchair

item #: 34860

sku : EIQ2N

brand : Sunrise Medical

\$1,230.00

List price: \$1,995.00

Leaves Warehouse In: 5-10 days

0 0

Quick Overview

- The best selling Ultralight Wheelchair in the world for over 25 years
- More features and options than any other folding wheelchair

Specifications

Weight Capacity	265 lbs.
Weight	29 lbs.
Seat Width	12" - 22"
Seat Depth	12" - 20"
Seat-to-Floor Height	13" - 22"
Overall Width	20" - 29"
Folded Width	13"
Frame Color	Black - Yellow
Front Wheel Size	3" - 8"
Rear Wheel Size	20" - 26"
Axle Type	Quick Release
Armrests	Desk-Length, Full-Length
Back Height	9" - 21"
Upholstery Type	Nylon
Upholstery Color	Black
Legrests	Swingaway, Elevating
Medicare HCPCS Code	K0005
Warranty on Frame	Lifetime Limited
Warranty on Wearable Parts	One Year Limited
Shipping Weight	50 lbs.
Box Length	33"
Box Width	13"
Box Height	38"